

# Spending too much time with your CBC analyzer?

## **Switch to Sight OLO®**

No routine cleaning. No wet reagents. No daily QC.

### Easy, maintenance-free lab-grade CBC results in minutes.



**19 parameters** with 5-part differential CBC



**No wet reagents** or liquid waste



Easy sample preparation



**2 drops of blood** from a finger prick or venous draw



QC is not a manufacturer requirement



**No operator calibration**: Installed factory-calibrated



Remote Servicing



Minimal training and **no memorization** required

Contact sales@sightdx.com or visit sightdx.com to learn more about easy, maintenance-free CBC.

#### No mess, no stress sample preparation.

Simplify blood sample preparation with OLO's single-use cartridges. Just two drops of blood (27µl) required.



# SIZhtolo®